


**43102** **NEW!**

One Port Drip Chamber Cover  
Fits 0.27 inch OD Tubing (6.8 mm)  
PVC



**41441**

One Port Drip Chamber Cover  
Fits 0.27 inch OD Tubing  
(6.8 mm)  
PVC - Rigid



**43100** **NEW!**

Two Port Drip Chamber Cover  
One Port Fits 0.27 inch OD Tubing,  
One Port Fits 0.22 inch OD Tubing  
(6.8 mm, 5.5 mm)  
PVC



**43203**

Two Port Drip Chamber Cover  
One Port Fits 0.27 inch OD Tubing,  
One Port Fits 0.22 inch OD Tubing  
(6.8 mm, 5.5 mm)  
PVC



**43204**

Soft Drip Chamber  
1.18 inch Top Diameter  
Fits 0.27 inch OD Tubing  
(30 mm; 6.8 mm)  
PVC



**41440**

Drip Chamber  
1.18 inch Top Diameter  
Fits 0.25 inch OD Tubing  
(30 mm; 6.3 mm)  
PVC



**NEW!**

**43101**

Drip Chamber  
1.18 inch Top Diameter  
Fits 0.27 inch - 0.25 inch OD Tubing  
(30 mm; 6.8 mm - 6.4 mm)  
PVC



Hand Pump



**43515**  
Drip Chamber, In-Line  
with Solution Filter  
15 Micron  
0.13 inch OD (3.3 mm)  
20 Drops/ml  
ABS, PA, SB



**43509**  
Drip Chamber, In-Line  
0.3 inch OD (3.3 mm)  
20 Drops/ml  
ABS, SB



**41445**  
Drip Chamber  
0.3 inch OD (3.3 mm)  
ABS



**41484**  
Hand Pump  
Fits 0.17 inch OD Tubing  
(4.3 mm)  
PP, PVC



Accessories

Chambers

Chamber Cover



**41443**  
One Port Drip Chamber Cover  
Fits 0.51 inch Drip Chamber (13 mm)  
0.16 inch ID x 0.25 inch OD  
(4 mm x 6.4 mm)  
ABS

Solution Filter



**X5004**  
Solution Filter  
15 Micron  
0.512 inch Diameter  
(13 mm)  
HDPE, Nylon



**23242**  
Solution Filter  
15 Micron  
0.52 inch Diameter  
(15.7 mm)  
ABS, PA



**X5001**  
Solution Filter  
15 Micron  
0.638 inch Diameter  
(16.2 mm)  
HDPE, Nylon



**X5002**  
Solution Filter  
25 Micron  
0.638 inch Diameter  
(16.2 mm)  
HDPE, Nylon



**43213**  
Solution Filter  
50 Micron  
0.638 inch Diameter  
(16.2 mm)  
ABS, Nylon

Infusion Filter



**X3001**  
Drip Chamber  
0.55 inch Top Diameter  
0.16 inch OD (14 mm; 4 mm)  
PVC - Soft



**23240**  
Drip Chamber  
0.63 inch Top Diameter  
0.16 inch ID (16 mm; 4 mm)  
PVC



**43206**  
Conical Filter  
Fits 0.63 inch Drip Chamber  
(16 mm)  
270 Micron  
PP



**41487**  
Filter  
Fits 6.7 inch Drip Chamber  
263 Micron  
2.44 inch Length  
(17 mm; 62 mm)  
Nylon, PP