




Guide Wire Dispenser



NEW!

11091

Fits 0.025 inch - 0.038 inch Guide Wires
(0.64 mm - 0.96 mm);
Accommodates 27 inch Guide Wire Length
(70 cm)
HDPE

Repositioning Sleeve



32200

0.119 inch ID;
31.5 inch
(3 mm; 80 cm)
ABS, LDPE, PE, PP

Guide Wire Dispensing Adapter

Female Luer Lock



11609

0.16 inch ID;
0.25 inch OD Tubing
(4.1 mm; 6.4 mm)
Acrylic