



J Straighteners



**32204**  
0.09 inch OD;  
Fits 0.025 - 0.03 inch Guide Wires  
(2.3 mm; 0.635 - 0.762 mm)  
PP



**32203**  
0.114 inch OD;  
Fits 0.032 - 0.038 inch Guide Wires  
(2.9 mm; 0.812 - 0.965 mm)  
PP



**11608**  
0.178 inch OD;  
Fits 0.039 inch Guide Wires  
(4.5 mm; 1 mm)  
PE, PP Blend




**11092** **NEW!**  
0.157 inch ID;  
Fits 0.039 inch Guide Wires  
(4 mm; 1 mm)  
MABS



**32201**  
0.114 inch OD;  
Fits 0.042 inch Guide Wires  
(2.9 mm; 1.06 mm)  
PP



**11630**  
0.178 inch OD;  
Fits 0.043 inch Guide Wires  
(4.5 mm; 1.09 mm)  
PE, PP Blend




**99911**  
0.169 - 0.153 inch ID;  
Fits 0.044 inch Guide Wires  
(4.3 - 3.9 mm; 1.11 mm)  
PE



**33030**  
0.114 inch OD;  
Fits 0.065 inch Guide Wires  
(2.9 mm; 1.65 mm)  
PP



**36203** **NEW!**  
0.039 inch - 0.044 inch ID  
(1 mm - 1.12 mm)  
PP



**36211**  
Base Plate for 36210 -  
Guide Wire Advancer  
ABS




**36210**  
0.061 inch ID (1.5 mm)  
LDPE


J Guide Wire Tip Protectors



**11090** **NEW!**  
0.82 inch (21 mm)  
PP



**32202**  
0.98 inch (25 mm)  
Dynaflex, Kraton - Blend



**32110**  
1.16 inch (29.5 mm)  
PE