





Female Luer Lock

Syringes



C3602
1.2 ml
PP, Santoprene



65732 **NEW!**
3 ml
PP, Silicone



C3603
5 ml
PP, Santoprene




C5005
50 ml
PP, Polyisoprene



11166
Non-Vented Male Luer Lock Cap
PP


Female Luer Lock Syringes can be attached to any Male Luer Lock Syringe, allowing the mixing of fluids between two syringes and eliminating the need for adapters.

Larger thru hole permits the passage of high viscosity fluids or pastes, as well as low viscosity materials.




C3608
1 ml
Polyisoprene
Rubber, PP

NEW!




C3610
5 ml
Polyisoprene
Rubber, PP

NEW!




C3612
20 ml
Polyisoprene
Rubber, PP

NEW!



C3609
2.5 ml
Polyisoprene
Rubber, PP


NEW!



C3611
10 ml
Polyisoprene
Rubber, PP

NEW!

Eccentric
Female
Luer Lock



C3613
60 ml
Polyisoprene
Rubber, PP

NEW!

Eccentric
Female
Luer Lock