



Male Luer Slip with Shroud



Shrouded Taper

80339

3 FR (0.039 inch, 1 mm)
20 Gauge
ABS, PP, Silicone



Shrouded Taper

80342

9 FR (0.118 inch, 3 mm)
11 Gauge
ABS, PP, Silicone

Rotating Male Hub Tuohy Borst with Sideport



Female Luer Lock with
Male Cap

80346

6.5 inch Length, 4 FR (0.053 inch ID, 1.35 mm) Gasket, 0.06 inch ID Cap (1.5 mm)
LDPE, MABS, PC, PP, PVC, Silicone