



Tuohy Borst Components

Tuohy Borst

It is only possible to combine Flat Cap, Extruded Gasket and Male Luer Lock Body components



NEW!

Part Number	Description	Size	Material
80416	Tuohy Borst Adapter Cap, Natural	0.165 inch Thru Hole (4.1 mm)	PE
80417	Tuohy Borst Adapter Cap, Blue	0.144 inch Thru Hole (3.7 mm)	ABS/Nylon
80418	Tuohy Borst Adapter Cap, Blue	0.078 inch Thru Hole (2 mm)	ABS/Nylon
80419	Tuohy Borst Adapter Cap, Blue	0.12 inch Thru Hole (3 mm)	ABS/Nylon




NEW!

Part Number	Description	Size	Material
80428	Tuohy Borst Adapter Gasket, Extruded Style	4 FR, 0.052 inch ID x 0.25 inch OD x 0.25 inch (1.5 mm x 6.3 mm x 6.3 mm)	Silicone
80429	Tuohy Borst Adapter Gasket, Extruded Style	9 FR, 0.113 inch ID x 0.25 inch OD x 0.156 inch (2.8 mm x 6.35 mm x 3.9 mm)	Silicone
80430	Tuohy Borst Adapter Gasket, Extruded Style	8 FR, 0.1 inch ID x 0.25 inch OD x 0.25 inch (2.5 mm x 6.35 mm x 6.35 mm)	Silicone
80431	Tuohy Borst Adapter Gasket, Extruded Style	6 FR, 0.081 inch ID x 0.25 inch OD x 0.25 inch (2 mm x 6.35 mm x 6.35 mm)	Silicone
80432	Tuohy Borst Adapter Gasket, Extruded Style	5.5 FR, 0.074 inch ID x 0.195 inch OD x 0.45 inch (1.8 mm x 4.9 mm x 11.3 mm)	Silicone
80433	Tuohy Borst Adapter Gasket, Extruded Style	4 FR, 0.055 inch ID x 0.19 inch OD x 0.406 inch (1.4 mm x 4.8 mm x 10 mm)	Silicone
80434	Tuohy Borst Adapter Gasket, Extruded Style	5 FR, 0.062 inch ID x 0.25 inch OD x 0.25 inch (1.59 mm x 6.35 mm x 6.35 mm)	Silicone
80435	Tuohy Borst Adapter Gasket, Extruded Style	3 FR, 0.04 inch ID x 0.25 inch OD x 0.25 inch (1 mm x 6.35 mm x 6.35 mm)	Silicone
80436	Tuohy Borst Adapter Gasket, Extruded Style	2 FR, 0.025 inch ID x 0.25 inch OD x 0.25 inch (0.64 mm x 6.35 mm x 6.35 mm)	Silicone



NEW!

Part Number	Description	Size	Material
80440	Tuohy Borst Adapter Body, Male Luer Lock with Side Port	0.125 inch ID (3.2 mm), 0.12 inch Thru Hole (3 mm)	PVC
80442	Tuohy Borst Adapter Body, Male Luer Lock with Side Port	0.125 inch ID (3.2 mm), 0.115 inch Thru Hole (2.9 mm)	PVC
80444	Tuohy Borst Adapter Body, Male Luer Lock 	0.12 inch Thru Hole (3 mm)	PVC
80445	Tuohy Borst Adapter Body, Male Luer Lock	0.12 inch Thru Hole (3 mm)	PVC



It is only possible to combine Female Luer Lock Cap, Molded Gasket and Snap Style Body components



NEW!

Part Number	Description	Size	Material
80421	Tuohy Borst Adapter Cap, Blue, Female Luer Lock	0.12 inch Thru Hole (3 mm)	ABS/Nylon
80422	Tuohy Borst Adapter Cap, Blue, Female Luer Lock with Probe Extension	0.112 inch Thru Hole (2.8 mm)	ABS/Nylon
80423	Tuohy Borst Adapter Cap, Clear, Female Luer Lock with Probe Extension	0.112 inch Thru Hole (2.8 mm)	PVC
80424	Tuohy Borst Adapter Cap, Red, Female Luer Lock	0.089 inch Thru Hole (2.3 mm)	ABS/Nylon
80425	Tuohy Borst Adapter Cap, Blue, Female Luer Lock	0.089 inch Thru Hole (2.3 mm)	ABS/Nylon
80426	Tuohy Borst Adapter Cap, Natural, Female Luer Lock	0.07 inch Thru Hole (1.78 mm)	ABS/Nylon



NEW!

Part Number	Description	Size	Material
80437	Tuohy Borst Adapter Gasket, Molded Style, Natural	6 FR, 0.084 inch ID (2.1 mm)	Kraton
80438	Tuohy Borst Adapter Gasket, Molded Style, Natural	4 FR, 0.057 inch ID (1.4 mm)	Kraton
80439	Tuohy Borst Adapter Gasket, Molded Style, Natural	3 FR, 0.039 inch ID (1 mm)	Kraton



NEW!

Part Number	Description	Size	Material
80446	Tuohy Borst Adapter Body, Snap Style	0.089 inch Thru Hole (2.3 mm)	PVC
80447	Tuohy Borst Adapter Body, Snap Style γ	0.055 inch Thru Hole (1.4 mm)	PVC