



Catheter Hub



41523
2.7 FR
(0.035 inch;
0.9 mm)
PC



51502
2.7 FR
(0.035 inch;
0.9 mm)
PC



41425
4 FR
(0.058 inch -
0.056 inch;
1.5 mm)
PC



41518
4 FR
(0.058 inch -
0.056 inch;
1.5 mm)
PC



41426
5 FR
(0.071 inch -
0.069 inch;
1.8 mm)
PC



41519
5 FR
(0.071 inch -
0.069 inch;
1.8 mm)
PC



41427
6 FR
(0.083 inch - 0.081 inch;
2.1 mm)
PC



41520
6 FR
(0.083 inch - 0.081 inch;
2.1 mm)
PC




41521
8 FR ID
(0.109 inch - 0.107 inch;
2.8 mm)
PC




41428
8 FR ID
(0.109 inch - 0.107 inch;
2.8 mm)
PC


Fixing Wings



41421
4.5 FR Thru Hole
(0.058 inch; 1.47 mm)
Pebax




65821
4.5 FR Thru Hole
(0.058 inch; 1.47 mm)
Kraton




65822
7 FR Thru Hole
(0.091 inch; 2.3 mm)
PVC



41423
7 FR Thru Hole
(0.091 inch; 2.3 mm)
Pebax



21500
12 FR Thru Hole
(0.156 inch; 4 mm)
HDPE



65199
14 FR Thru Hole
(0.180 inch; 4.5 mm)
PVC/PU Blend



65198
15 FR Thru Hole
(0.195 inch; 5 mm)
PVC/PU Blend