



Chamber Cover



41443
One Port Drip Chamber Cover
Fits 0.51 inch Drip Chamber (13 mm)
0.16 inch ID x 0.25 inch OD
(4 mm x 6.4 mm)
ABS

Solution Filter



X5004
Solution Filter
15 Micron
0.512 inch Diameter
(13 mm)
HDPE, Nylon



23242
Solution Filter
15 Micron
0.52 inch Diameter
(15.7 mm)
ABS, PA



X5001
Solution Filter
15 Micron
0.638 inch Diameter
(16.2 mm)
HDPE, Nylon



X5002
Solution Filter
25 Micron
0.638 inch Diameter
(16.2 mm)
HDPE, Nylon



43213
Solution Filter
50 Micron
0.638 inch Diameter
(16.2 mm)
ABS, Nylon

Infusion Filter



X3001
Drip Chamber
0.55 inch Top Diameter
0.16 inch OD (14 mm; 4 mm)
PVC - Soft



23240
Drip Chamber
0.63 inch Top Diameter
0.16 inch ID (16 mm; 4 mm)
PVC



43206
Conical Filter
Fits 0.63 inch Drip Chamber
(16 mm)
270 Micron
PP



41487
Filter
Fits 6.7 inch Drip Chamber
263 Micron
2.44 inch Length
(17 mm; 62 mm)
Nylon, PP