




Shown Assembled

Sanitary Series Couplings feature an ergonomic thumb latch which is easy to operate, even by gloved hands. The 3/4 inch and 1 inch sanitary couplings can be installed with standard gaskets and clamps. This series is also compatible with the MPC and MPX series, which allows for easy conversion to industry standard connections or single-use systems.



A = Height/Diameter  
B = Total Length (Straight Thru/Shut Off)  
E = Outside Diameter

### Coupling Bodies




In-Line  
PSU

In-Line Sanitary	Size	A x B x E
<b>MPC3301239</b>	3/4 inch Mini (19.1 mm)	0.98 inch x 1.4 inch x 1 inch (24.9 mm x 35.6 mm x 25.4 mm)



In-Line  
PSU

In-Line Sanitary	Size	A x B x E
<b>MPC3301639</b>	1 inch Maxi (25.4 mm)	1.5 inch x 1.4 inch x 1.5 inch (38.1 mm x 35.6 mm x 38.1 mm)



In-Line  
PSU

In-Line Sanitary	Size	A x B x E
<b>MPX3301239</b>	3/4 inch Mini (19.1 mm)	1.28 inch x 1.7 inch x 1 inch (32.5 mm x 43.2 mm x 25.4 mm)