



Barbed & Stepped Connectors

Connectors



21060
Both Ports 0.33 inch -
0.46 inch OD
(8.4 mm - 11.7 mm)
PVC



21070
0.33 inch - 0.47 inch OD x
0.38 inch - 0.43 inch OD
(8.4 mm - 11.9 mm x
9.7 mm - 10.9 mm)
PVC



21075
0.51 inch - 0.59 inch OD x
0.63 inch OD
(13 mm - 15 mm x 16 mm)
PVC



21068
Fits 0.351 inch - 0.522 inch ID x
0.363 inch ID Tubing
(8.9 mm - 13.3 mm x 9.2 mm)
PS



21177
Fits 0.212 inch - 0.358 inch ID Tubing x
0.315 inch ID
(5.4 mm - 9.1mm x 8 mm)
PVC - Soft



11179
0.39 inch - 0.55 inch OD
(10.1 mm - 13.8 mm)
ABS



11181
0.28 inch - 0.43 inch OD
(7.1 mm - 10.9 mm)
ABS



94827
Vented
PP



94627
Non-Vented
Fits 0.259 inch - 0.394 inch ID Tubing x
0.197 inch OD
(6.6 mm - 10 mm x 5 mm)
PC



94817
Non-Vented
PP