



28224
Epidural Filter with Female Luer Lock Inlet
to Male Luer Lock Outlet
0.2 Micron
Acrylic, Polyester - Non-Woven, Supor PES

Front View



Back View



28225
Neonatal/Pediatric IV Filter
0.2 Micron
0.138 inch OD (3.5 mm)
Acrylic, PES, PTFE



28222
Micro-Infusion IV Filter
0.2 Micron
0.08 inch - 0.05 inch ID (2 mm - 1.3 mm)
Acrylic, PES



12114
Hydrophilic In-Line IV Filter
0.22 Micron
0.09 inch OD (2.3 mm)
Acrylic, PES, PVDF