



11679
Insufflation Filter with Hose Barb
0.1 Micron
50 mm Diameter, 0.34 inch Hose Barb (8.6 mm)
K-Resin, PTFE



28205
Hydrocheck Gas Line Filter with Graduated Hose Barb
0.1 Micron
57 mm Diameter, 1/4 inch Graduated Hose Barb
(0.25 inch; 6.4 mm)
Acrylic, PTFE



28204
Hydrocheck Gas Line Filter with Graduated Hose Barb
0.1 Micron
57 mm Diameter, 1/4 inch Graduated Hose Barb
(0.25 inch; 6.4 mm)
Acrylic, Microglass



28241 **NEW!**
Economy Gas Filter (PTFE) with Graduated Hose Barb
0.3 Micron
62 mm Diameter, 1/4 inch - 1/2 inch Graduated Hose Barb
(0.25 inch - 0.5 inch; 6.4 mm - 12.7 mm)
Acrylic, PTFE



28242 **NEW!**
Economy Gas Filter (Microglass) with Graduated Hose Barb
0.3 Micron
62 mm Diameter, 1/4 inch - 1/2 inch Graduated Hose Barb
(0.25 inch - 0.5 inch; 6.4 mm - 12.7 mm)
Acrylic, Microglass