



In-Line

Filters



28233
Hydrophilic Filter Hub with
Female Luer Lock to
Male Luer Lock
5 Micron
ABS, Nylon



28232
Hydrophilic Filter Hub with
Female Luer Lock to
Male Luer Slip
5 Micron
ABS, Nylon




12113
Hydrophilic In-Line IV Filter,
Female Luer Lock Inlet
0.22 Micron
0.13 inch OD (3.3 mm)
COPE, PETG, PES, PDVF



12115
Hydrophilic In-Line IV Filter,
Female Luer Lock Inlet
1.2 Micron
0.13 inch OD (3.3 mm)
COPE, PETG, PES, PDVF


Flow Direction →



Male Luer Slip


NEW!

28236
Hydrophilic Filter with Female Luer Lock Inlet
to Male Luer Slip Outlet
0.2 Micron
13 mm
CN-CA, PP



Male Luer Slip


28235
Hydrophilic Filter with Female Luer Lock Inlet
to Male Luer Slip Outlet
0.2 Micron
20 mm
Nylon, PP



Male Luer Slip

Female Luer Slip

28215
Hydrophobic Filter with Male Luer Slip Inlet
to Female Luer Slip Outlet
0.2 Micron
25 mm
Acrylic, PTFE




Male Luer Slip


28217
Hydrophilic Filter with Female Luer Lock Inlet
to Male Luer Slip Outlet
0.2 Micron
25 mm
Acrylic, Supor




Flow Direction →




28219
Hydrophobic Filter with Female Luer Slip Inlet
to Male Luer Slip Outlet
0.2 Micron
17 mm
Acrylic, PTFE




28213
Hydrophobic Filter with Female Luer Lock Inlet
to Male Luer Lock Outlet
0.2 Micron
17 mm
Acrylic, PTFE




28203
Hydrophobic Filter with Female Luer Lock Inlet
to Male Slip Outlet
0.2 Micron
25 mm
PVC, PTFE



28216
Hydrophilic Filter with Female Luer Lock Inlet
to Male Luer Lock Outlet
0.2 Micron
25 mm
Acrylic, Supor



28234 **NEW!**
Hydrophobic Filter with Female Luer Lock Inlet
to Male Luer Slip Outlet
0.2 Micron
25 mm
PP, PVDF



28211
Hydrophobic Filter with Female Luer Lock Inlet
to Male Luer Lock Outlet
0.45 Micron
17 mm
Acrylic, PTFE