



Flow Direction →



**28220**  
Adult IV Filter, Self Priming  
0.2 Micron  
0.12 inch - 0.125 inch OD (3 mm - 3.2 mm)  
Acrylic, PES, PTFE



**28218**  
Adult IV Filter  
0.2 Micron  
0.12 inch - 0.125 inch OD (3 mm - 3.2 mm)  
Acrylic, Supor



**28230**  
Adult IV Filter, Positively Charged  
0.2 Micron  
0.12 inch OD (3 mm)  
Acrylic, PES, PTFE



**28229**  
Adult IV Filter, Self Priming  
0.2 Micron  
0.15 inch ID (3.8 mm)  
Acrylic, PES, PTFE



**28221**  
Adult IV Filter, Self Priming  
1.2 Micron  
0.12 inch - 0.135 inch OD (3 mm - 3.4 mm)  
Acrylic, PES, PTFE



**13126**  
Adult IV Filter, Self Priming  
1.2 Micron  
0.12 inch OD (3 mm)  
Acrylic, PES, PTFE



**28224**  
Epidural Filter with Female Luer Lock Inlet  
to Male Luer Lock Outlet  
0.2 Micron  
Acrylic, Polyester - Non-Woven, Supor PES

Front View



Back View



**28225**  
Neonatal/Pediatric IV Filter  
0.2 Micron  
0.138 inch OD (3.5 mm)  
Acrylic, PES, PTFE



**28222**  
Micro-Infusion IV Filter  
0.2 Micron  
0.08 inch - 0.05 inch ID (2 mm - 1.3 mm)  
Acrylic, PES



**12114**  
Hydrophilic In-Line IV Filter  
0.22 Micron  
0.09 inch OD (2.3 mm)  
Acrylic, PES, PVDF