


Female Luer Slip Female Luer Slip

32101


Transducer Protector
without Adapter
Fits 0.22 inch OD Tubing (5.6 mm)
Hydrophobic, PVC



Female Luer Slip

32103


Transducer Protector
with Female Luer Lock and Tube Site
Fits 0.22 inch OD Tubing (5.6 mm)
Hydrophobic, PVC



Female Luer Slip Female Luer Slip

32106


Transducer Protector
with Double Tube Site
Fits 0.22 inch OD Tubing (5.6 mm)
Hydrophobic, PVC



Female Luer Slip Female Luer Slip

11500


Transducer Protector, Blue
Fits 0.22 inch OD Tubing (5.6 mm)
Hydrophobic, PVC



Female Luer Slip Female Luer Slip

11501


Transducer Protector, Red
Fits 0.22 inch OD Tubing (5.6 mm)
Hydrophobic, PVC



Male Luer Slip Male Luer Slip


32105

Transducer Protector
Hydrophobic, PVC



32111


Transducer Protector
with Female Luer Lock Inlet
and Male Luer Lock Outlet
Hydrophobic, PVC



Male Luer Slip

32107


Transducer Protector
with Male Luer Slip
and Male Luer Lock
Hydrophobic, PVC



32104

Transducer Protector
with Female and Male Luer Lock,
Clear/White
Hydrophobic, PVC

Snap-Off Male Cap for
Female Luer Lock




Male Luer Slip

32108

Transducer Protector
with Male Luer Slip,
Features Snap-Off Cap on Strap
Hydrophobic, PVC


Snap-Off Male Cap for
Female Luer Lock



32109

Transducer Protector
with Male Luer Lock,
Features Snap-Off Cap on Strap
Hydrophobic, PVC

Snap-Off Male Cap for
Female Luer Lock



32102

Transducer Protector
with Snap off Male Cap for
Female Luer Lock, Tubing Port
Fits 0.22 inch OD Tubing (5.6 mm)
Hydrophobic, PVC