



Non-Latex Injection Sites

Injection Sites

Male Luer Injection Sites



800211
PC, Polyisoprene



Spin Lock
55906
ABS, Polyisoprene
Pulls Down



80181 **NEW!**
ABS, Polyisoprene



91015
ABS, Isoprene



80028 **NEW!**
PC, Polyisoprene



80029
PC, Polyisoprene



80034
PC, Polyisoprene

T Injection Sites



80016
0.094 inch ID (2.4 mm)
ABS, Isoprene, SB, PVC



80204
0.138 inch ID (3.5 mm)
ABS, Polyisoprene



Y Injection Sites



80203

0.079 inch ID (2 mm) Both Ports
ABS, Polyisoprene



Male Luer Slip

12103

0.086 inch ID (2.2 mm) Both Ports,
0.138 inch OD (3.5 mm) Side Port
ABS, Polyisoprene




80012

0.156 inch ID (4.0 mm) Both Ports
ABS, Polyisoprene



12102

0.157 inch ID (4 mm) Both Ports
ABS, Polyisoprene



80005

0.161 inch ID (4.1 mm) Both Ports
ABS, Polyisoprene, PVC



12106

Fits 0.268 inch ID Tubing
(6.8 mm)
ABS, Polyisoprene



84033

PC, Polyisoprene



84037

0.080 inch ID x 0.125 inch OD
(2 mm x 3.2 mm) Both Ports
Polyisoprene, PC



84038

0.116 inch ID x 0.160 inch OD;
0.080 inch ID x 0.125 inch OD
(2.9 mm x 4.1 mm; 2 mm x 3.2 mm)
Polyisoprene, PC

Non-Latex Injection Sites



Injection Sites



49009

0.3 inch - 0.276 inch OD (5.6 mm - 7 mm)
Bromobutyl, PC, PP



49005

0.189 inch ID x 0.287 inch OD (4.8 mm x 6.8 mm)
Al, Latex Rubber, PC, PP



80055

0.251 inch OD x 0.236 inch OD
(6.4 mm x 6 mm)
Polyisoprene



11636

Fits 0.268 inch ID Port (6.8 mm)
PC, Polyisoprene


Break-Off Tip

Barbed



11456

Fits 1/16 inch ID Tubing
(1.6 mm)
Isoplast, Polyisoprene, PP



11194

Fits 1/16 inch ID Tubing
(1.6 mm)
Isoplast, Polyisoprene, PP

Pre-Slit



80137

Fits 0.236 inch ID Tubing (6 mm)
PC, PP, Polyisoprene



80138

Fits 0.236 inch ID Tubing (6 mm)
PC, PP, Polyisoprene



80139

Fits 0.236 inch ID Tubing (6 mm)
PC, PP, Polyisoprene