



29205
4 ml Capacity, 3 ml Draw, 0.098 inch Stem Diameter (2.5 mm)
3.3 inch (85 mm)
LDPE



29215
1.7 ml Capacity, 0.9 ml Draw, 0.189 inch Stem Diameter (4.8 mm)
3.4 inch (87 mm)
LDPE




29201
1 ml Capacity, 0.8 ml Draw, 0.06 inch Stem Diameter (1.3 mm)
4.3 inch (109 mm)
PE



29236
4 ml Capacity, 1.9 ml Draw, 0.26 inch Stem Diameter (6.6 mm)
5 inch (129 mm)
LDPE



29230
2 ml Capacity, 1 ml Draw, 0.15 inch Stem Diameter (3.8 mm)
5.3 inch (137 mm)
LDPE



29190
7 ml Capacity, 2 ml Draw, 0.32 inch Stem Diameter (8.15 mm)
5.3 inch (137 mm)
LDPE



29210
3.9 ml Capacity, 2.2 ml Draw, 0.185 inch Stem Diameter (4.7 mm)
5.5 inch (140 mm)
LDPE




29235
4.6 ml Capacity, 3.3 ml Draw, 0.106 inch Stem Diameter (2.7 mm)
6.1 inch (157 mm)
LDPE



29200
5.8 ml Capacity, 3.4 ml Draw, 0.181 inch Stem Diameter (4.6 mm)
6.1 inch (157 mm)
LDPE



29220
7.5 ml Capacity, 3.2 ml Draw, 0.311 inch Stem Diameter (7.9 mm)
6.1 inch (157 mm)
LDPE



29191
15 ml Capacity, 5 ml Draw, 0.39 inch Stem Diameter (9.9 mm)
7.5 inch (190 mm)
LDPE