



Double Dead-Enders

Luers

Each cap includes both male and female ends.
Decreases the type and number of caps in a set.
Allows for more efficient compact packaging.

Vented



11455
PC



12094
PP



17659
HDPE



17662
HDPE



65835
PC



65312
PC



65313
PC



65314
PC



65813
PC




65814
PC



65815
PC

Non-Vented Recessed Male Luer Locks

Ideal for Luer Activated Valves. Guards against touch contamination.



65800
PP



65801
PP



65802
PP



65807
PP



65833
PP



65834
PP



Each cap includes both male and female ends.
Decreases the type and number of caps in a set.
Allows for more efficient compact packaging.

Non-Vented



NEW!

11086
HDPE



12092
PP



12093
PP



13815
PC



17660
HDPE



17661
HDPE



65804
PC



65811
PC



65812
PC



65826
PE