



Male Luer Slips

Luers

Barbed



NEW!

11734
Fits 1/32 inch
ID Tubing
(0.031 inch;
0.8 mm)
PC



89346
Fits 1/16 inch
ID Tubing
(0.062 inch;
1.6 mm)
ABS



89352
Fits 1/16 inch
ID Tubing
(0.062 inch;
1.6 mm)
PC



11660
Fits 1/16 inch
ID Tubing
(0.062 inch;
1.6 mm)
ABS



11666
Fits 1/16 inch
ID Tubing
(0.062 inch;
1.6 mm)
PP



11663
Fits 1/16 inch
ID Tubing
(0.062 inch;
1.6 mm)
PC



11662
Fits 3/32 inch
ID Tubing
(0.093 inch;
2.4 mm)
ABS



11665
Fits 3/32 inch
ID Tubing
(0.093 inch;
2.4 mm)
PC



11668
Fits 3/32 inch
ID Tubing
(0.093 inch;
2.4 mm)
PP



89347
Fits 3/32 inch
ID Tubing
(0.093 inch;
2.4 mm)
ABS



89353
Fits 3/32 inch
ID Tubing
(0.093 inch;
2.4 mm)
PC



11661
Fits 1/8 inch
ID Tubing
(0.125 inch;
3.2 mm)
ABS



11664
Fits 1/8 inch
ID Tubing
(0.125 inch;
3.2 mm)
PC



11667
Fits 1/8 inch
ID Tubing
(0.125 inch;
3.2 mm)
PP



89348
Fits 1/8 inch
ID Tubing
(0.125 inch;
3.2 mm)
ABS



89354
Fits 1/8 inch
ID Tubing
(0.125 inch;
3.2 mm)
PC



89349
Fits 5/32 inch
ID Tubing
(0.156 inch;
4 mm)
ABS



89355
Fits 5/32 inch
ID Tubing
(0.156 inch;
4 mm)
PC



89350
Fits 3/16 inch
ID Tubing
(0.187 inch;
4.7 mm)
ABS



89356
Fits 3/16 inch
ID Tubing
(0.187 inch;
4.7 mm)
PC



89351
Fits 1/4 inch
ID Tubing
(0.25 inch;
6.4 mm)
ABS



89357
Fits 1/4 inch
ID Tubing
(0.25 inch;
6.4 mm)
PC



11111
0.3 inch OD
Barbed End
(7.6 mm)
PVC



X7001
0.071 inch ID
x 0.157 inch -
0.181 inch OD
(1.8 mm x 4 mm
- 4.6 mm)
PVC



11117
0.105 inch - 0.071 inch ID
(2.66 mm - 1.8 mm)
PVC



20017
0.125 inch - 0.129 inch OD
(3.2 mm - 3.3 mm)
ABS



11167
0.145 inch ID
(3.7 mm)
PVC



71288
0.145 inch -
0.134 inch ID x
0.225 inch OD
(3.7 mm - 3.4 mm x
5.7 mm)
Acrylic



20013
0.157 inch - 0.149 inch ID x
0.25 inch - 0.254 inch OD
(4 mm - 3.8 mm x 6.4 mm -
6.5 mm)
ABS



11118
0.161 inch ID
(4.1 mm)
ABS



12108
0.161 inch ID
(4.1 mm)
ABS



11102
0.21 inch ID
(5.3 mm)
PVC

Male Slip to Male Slip



20036
ABS



11112
PVC