



Metal Luer Connectors

Luers

Male Luer to ID Tubing Connectors



99154
Fits 1/16 inch -
3/32 inch ID Tubing
(0.062 inch - 0.093 inch;
1.6 mm - 2.4 mm)
Nickel Plated Brass



99156
Fits 1/16 inch -
3/32 inch ID Tubing (0.062 inch -
0.093 inch;
1.6 mm - 2.4 mm)
Nickel Plated Brass



99152
Fits 1/8 inch -
3/16 inch ID Tubing
(0.125 inch - 0.187 inch;
3.2 mm - 4.7 mm)
Nickel Plated Brass



99153
Fits 1/8 inch -
3/16 inch ID Tubing
(0.125 inch - 0.187 inch;
3.2 mm - 4.7 mm)
Nickel Plated Brass



99155
Fits 1/8 inch - 3/16 inch ID Tubing
(0.125 inch - 0.187 inch;
3.2 mm - 4.7 mm)
Nickel Plated Brass



99150
Fits 1/4 inch ID Tubing
(0.25 inch; 6.4 mm)
Nickel Plated Brass



99151
Fits 1/4 inch ID Tubing
(0.25 inch; 6.4 mm)
Nickel Plated Brass



99160
Fits 1/4 inch -
5/16 inch ID Tubing
(0.25 inch - 0.312 inch;
6.4 mm - 7.9 mm)
Nickel Plated Brass

Female Luer Lock to ID Tubing Connectors



99158
Fits 1/16 inch - 3/32 inch ID Tubing
(0.062 inch - 0.093 inch;
1.6 mm - 2.4 mm)
Nickel Plated Brass



99157
Fits 1/8 inch - 3/16 inch ID Tubing
(0.125 inch - 0.187 inch;
3.2 mm - 4.7 mm)
Nickel Plated Brass



99159
Fits 3/16 inch - 1/4 inch ID Tubing
(0.187 inch - 0.25 inch;
4.7 mm - 6.4 mm)
Nickel Plated Brass



99161
Female Luer Lock
to Male Luer Slip
Nickel Plated Brass



99162
Female Luer Lock
to Male Luer Slip
Nickel Plated Brass



99163
Female Luer Lock
to Male Record
Nickel Plated Brass



99164
Female Record
to Male Luer Slip
Nickel Plated Brass

Connectors are also available in Stainless Steel and customized configurations