



Y Connectors with Luers

Luers



**80389**  
Top Port Female Luer Lock;  
Side Arm & Bottom Ports  
0.079 inch ID x 0.157 inch OD  
(2 mm x 4 mm)  
ABS




**83015**  
Side Arm & Top Ports  
Female Luer Locks;  
Bottom Port 0.118 inch ID  
(3 mm)  
PC



**83010**  
Side Arm & Top Ports  
Female Luer Locks;  
Bottom Port 0.125 inch ID  
(3.2 mm)  
ABS



**71622**  
Side Arm Female Luer Lock;  
Top & Bottom Ports  
0.146 inch - 0.135 inch ID x  
0.254 inch OD  
(3.7 mm - 3.4 mm x 6.5 mm)  
Acrylic



**65900**  
Side Arm  
Female Luer Lock;  
Top & Bottom Ports  
0.157 inch ID (4 mm)  
PVC



**11658**  
Side Arm & Top Ports  
Female Luer Locks;  
Bottom Port  
0.16 inch ID (4.1 mm)  
PVC



**80056**  
Side Arm & Top Ports  
Female Luer Locks;  
Bottom Port Male Slip  
ABS



**80413** **NEW!**  
Top Ports Male Slips;  
Bottom Port  
Female Luer Lock  
ABS



**84044**  
Side Arm & Top Ports  
Female Luer Locks;  
Bottom Port Male Luer Spin Lock  
ABS, PC




**84049** **NEW!**  
Top Ports  
Female Luer Locks;  
Bottom Port  
Male Luer Lock  
ABS, PP



**80155**  
Side Arm & Top Ports  
Female Slips;  
Bottom Port Male Slip  
PVC



**80160**  
Side Arm & Top Ports  
Female Luer Locks,  
Bottom Port  
0.225 inch OD (5.7 mm)  
PVC




**83017**  
Side Arm & Top Ports  
Female Luer Locks;  
Bottom Port Male Slip  
PC



**11454**  
Side Arm & Top Ports  
Male Luer Spin Locks;  
Bottom Port Female Luer Connector  
ABS, PC




**83009**  
Side Arm & Top Ports  
0.125 inch ID (3.2 mm);  
Bottom Port Male Luer Spin Lock  
ABS



**83016**  
Side Arm & Top Ports  
Female Luer Locks;  
Bottom Port Male Luer Spin Lock  
PC



**65925**  
Top Ports Female Luer Locks;  
Bottom Port Male Luer Lock  
PC, PVC



**80157**  
Top Ports Female Luer Locks,  
Bottom Port Male Luer Slip  
PVC