



Part Number	Description	Size	Material
<b>31414</b>	Needle Hub, Round with Female Luer Lock	14 G (0.079 inch, 2 mm)	PC
<b>31415</b>	Needle Hub, Round with Female Luer Lock	15 G (0.071 inch, 1.8 mm)	PC
<b>31416</b>	Needle Hub, Round with Female Luer Lock	16 G (0.063 inch, 1.6 mm)	PC
<b>31417</b>	Needle Hub, Round with Female Luer Lock, Universal	Fits Bushings 31217 - 31227	PC



Part Number	Description	Size	Material
<b>31514</b>	Needle Hub, Anatomic with Female Luer Lock	14 G (0.079 inch, 2 mm)	PC
<b>31515</b>	Needle Hub, Anatomic with Female Luer Lock	15 G (0.071 inch, 1.8 mm)	PC
<b>31516</b>	Needle Hub, Anatomic with Female Luer Lock	16 G (0.063 inch, 1.6 mm)	PC
<b>31517</b>	Needle Hub, Anatomic with Female Luer Lock, Universal	Fits Bushings 31217 - 31227	PC



Part Number	Description	Size	Material
<b>31714</b>	Needle Hub, Anatomic Octagonal with Female Luer Lock	14 G (0.079 inch, 2 mm)	PETG
<b>31715</b>	Needle Hub, Anatomic Octagonal with Female Luer Lock	15 G (0.071 inch, 1.8 mm)	PETG
<b>31716</b>	Needle Hub, Anatomic Octagonal with Female Luer Lock	16 G (0.063 inch, 1.6 mm)	PETG
<b>31717</b>	Needle Hub, Anatomic Octagonal with Female Luer Lock, Universal	Fits Bushings 31217 - 31227	PETG



NEW!

Part Number	Description	Size	Material
31814	Needle Hub with Wings and Female Luer Lock	14 G (0.079 inch, 2 mm)	PC
31815	Needle Hub with Wings and Female Luer Lock	15 G (0.071 inch, 1.8 mm)	PC
31816	Needle Hub with Wings and Female Luer Lock	16 G (0.063 inch, 1.6 mm)	PC
31817	Needle Hub with Wings and Female Luer Lock, Universal	Fits Bushings 31217 - 31227	PC



NEW!

Part Number	Description	Size	Material
31914	Needle Hub with Wings, Female Luer Lock to Male Luer Slip	14 G (0.079 inch, 2 mm)	PC
31915	Needle Hub with Wings, Female Luer Lock to Male Luer Slip	15 G (0.071 inch, 1.8 mm)	PC
31916	Needle Hub with Wings, Female Luer Lock to Male Luer Slip	16 G (0.063 inch, 1.6 mm)	PC
31917	Needle Hub with Wings, Female Luer Lock to Male Luer Slip, Universal	Fits Bushings 31217 - 31227	PC



NEW!

Part Number	Description	Material
31918	Spin Lock Collar for Needle Hub 31814 - 31817	PC



NEW!

Part Number	Description	Size	Material
36314	Needle Hub, Rectangular with Female Luer Lock	14 G (0.079 inch, 2 mm)	PC
36315	Needle Hub, Rectangular with Female Luer Lock	15 G (0.071 inch, 1.8 mm)	PC
36316	Needle Hub, Rectangular with Female Luer Lock	16 G (0.063 inch, 1.6 mm)	PC
36317	Needle Hub, Rectangular with Female Luer Lock, Universal	Fits Bushings 31217 - 31227	PC



Part Number	Description	Tubing Size	Material
31014	Cap, Dark Green	Fits 14 G Round Needle Hub	ABS
31015	Cap, Dark Blue	Fits 15 G Round Needle Hub	ABS
31016	Cap, White	Fits 16 G Round Needle Hub	ABS
31017	Cap, Red	Fits 17 G Round Needle Hub	ABS
31018	Cap, Pink	Fits 18 G Round Needle Hub	ABS
31019	Cap, Beige	Fits 19 G Round Needle Hub	ABS
31020	Cap, Yellow	Fits 20 G Round Needle Hub	ABS
31021	Cap, Light Green	Fits 21 G Round Needle Hub	ABS
31022	Cap, Black	Fits 22 G Round Needle Hub	ABS
31023	Cap, Light Blue	Fits 23 G Round Needle Hub	ABS
31024	Cap, Violet	Fits 24 G Round Needle Hub	ABS
31025	Cap, Orange	Fits 25 G Round Needle Hub	ABS
31026	Cap, Brown	Fits 26 G Round Needle Hub	ABS
31027	Cap, Grey	Fits 27 G Round Needle Hub	ABS



Part Number	Description	Tubing Size	Material
31614	Cap, Dark Green	Fits 14 G Anatomic Octagonal Needle Hub	PC
31615	Cap, Dark Blue	Fits 15 G Anatomic Octagonal Needle Hub	PC
31616	Cap, White	Fits 16 G Anatomic Octagonal Needle Hub	PC
31617	Cap, Red	Fits 17 G Anatomic Octagonal Needle Hub	PC
31618	Cap, Pink	Fits 18 G Anatomic Octagonal Needle Hub	PC
31619	Cap, Beige	Fits 19 G Anatomic Octagonal Needle Hub	PC
31620	Cap, Yellow	Fits 20 G Anatomic Octagonal Needle Hub	PC
31621	Cap, Light Green	Fits 21 G Anatomic Octagonal Needle Hub	PC
31622	Cap, Black	Fits 22 G Anatomic Octagonal Needle Hub	PC
31623	Cap, Light Blue	Fits 23 G Anatomic Octagonal Needle Hub	PC
31624	Cap, Violet	Fits 24 G Anatomic Octagonal Needle Hub	PC
31625	Cap, Orange	Fits 25 G Anatomic Octagonal Needle Hub	PC
31626	Cap, Brown	Fits 26 G Anatomic Octagonal Needle Hub	PC
31627	Cap, Grey	Fits 27 G Anatomic Octagonal Needle Hub	PC



NEW!

Part Number	Description	Tubing Size	Material
35314	Cap, Dark Green	Fits 14 G Rectangular Needle Hub	ABS
35315	Cap, Dark Blue	Fits 15 G Rectangular Needle Hub	ABS
35316	Cap, White	Fits 16 G Rectangular Needle Hub	ABS
35317	Cap, Red	Fits 17 G Rectangular Needle Hub	ABS
35318	Cap, Pink	Fits 18 G Rectangular Needle Hub	ABS
35319	Cap, Beige	Fits 19 G Rectangular Needle Hub	ABS
35320	Cap, Yellow	Fits 20 G Rectangular Needle Hub	ABS
35321	Cap, Light Green	Fits 21 G Rectangular Needle Hub	ABS
35322	Cap, Black	Fits 22 G Rectangular Needle Hub	ABS
35323	Cap, Light Blue	Fits 23 G Rectangular Needle Hub	ABS
35324	Cap, Violet	Fits 24 G Rectangular Needle Hub	ABS
35325	Cap, Orange	Fits 25 G Rectangular Needle Hub	ABS
35326	Cap, Brown	Fits 26 G Rectangular Needle Hub	ABS
35327	Cap, Grey	Fits 27 G Rectangular Needle Hub	ABS



NEW!

Part Number	Description	Size	Material
31217	Bushing to obtain 17 G	(0.058 inch, 1.47 mm)	PETG
31218	Bushing to obtain 18 G	(0.05 inch, 1.27 mm)	PETG
31219	Bushing to obtain 19 G	(0.043 inch, 1.1 mm)	PETG
31220	Bushing to obtain 20 G	(0.035 inch, 0.9 mm)	PETG
31221	Bushing to obtain 21 G	(0.031 inch, 0.8 mm)	PETG
31222	Bushing to obtain 22 G	(0.027 inch, 0.7 mm)	PETG
31223	Bushing to obtain 23 G	(0.023 inch, 0.6 mm)	PETG
31224	Bushing to obtain 24 G	(0.021 inch, 0.55 mm)	PETG
31225	Bushing to obtain 25 G	(0.02 inch, 0.5 mm)	PETG
31226	Bushing to obtain 26 G	(0.018 inch, 0.45 mm)	PETG
31227	Bushing to obtain 27 G	(0.016 inch, 0.4 mm)	PETG