

High-Flow Connectors



Non-LuERS




31340 **NEW!**
0.27 inch ID x 0.35 inch OD
(6.8 mm x 9 mm)
PVC



91001
0.26 inch - 0.252 inch ID x
0.354 inch - 0.339 inch OD
(6.8 mm - 6.5 mm x 9 mm - 8.6 mm)
PVC




49011 **NEW!**
Fits 0.27 inch OD Tubing
(6.8 mm)
PVC




91000
0.264 inch - 0.256 inch ID x
0.354 inch OD
(6.8 mm - 6.5 mm x 9 mm)
PVC



31318
Fits 0.268 inch OD
Tubing (6.8 mm)
PVC - Soft




31319
Fits 0.268 inch OD
Tubing (6.8 mm)
PVC - Soft



31320
Fits 0.268 inch OD
Tubing (6.8 mm)
PVC - Rigid



31321
Fits 0.268 inch OD
Tubing (6.8 mm)
PVC - Rigid



31322
Fits 0.268 inch OD
Tubing (6.8 mm)
PVC - Rigid

Hansen Connectors



28115
Fits 0.268 inch OD Tubing
(6.8 mm)
PVC



28116
Fits 0.268 inch OD Tubing
(6.8 mm)
PVC



28117
Fits 0.268 inch OD Tubing
(6.8 mm)
PVC