




Tear-off/Break-off Ports



Break-Off Tip

**89310**


0.213 inch ID x 0.28 inch OD  
(5.4 mm x 7.1 mm)  
PVC



Break-Off Tip

**63004**

Fits 0.252 inch ID Tubing x 0.276 inch OD  
(6.4 mm x 7 mm)  
PVC




Break-Off Tip

**NEW!**

**65840**

Bag Port, Twist-Off with Filter  
3 Micron  
0.197 inch ID x 0.276 inch OD  
(5 mm x 7 mm)  
Acrylic Copolymer, PVC



Break-Off Tip

**65823**

0.280 inch OD (7.1 mm)  
PVC



Break-Off Tip

**63017**

PVC



Break-Off Tip

**12105**

0.268 inch ID (6.8 mm)  
PVC



Break-Off Tip

**94284**

0.102 ID x 0.193 inch - 0.209 inch OD  
(2.6 mm x 4.9 mm - 5.3 mm)  
PC

Female Luer Lock  
with Cap Included



Break-Off Tip

**11503**

0.256 inch - 0.264 inch OD Tapered  
(6.5 mm - 6.7 mm)  
PC



Break-Off Tip

**94282**

Fits 1/4 inch ID Tubing (0.25 inch; 6.4 mm);  
Adapter: 0.268 inch - 0.283 inch OD  
(6.8 mm - 7.2 mm )  
PC