




Tracheal Tube Components

Respiratory



**53501**  
0.6 inch x 0.1 inch  
(15 mm x 2.5 mm)  
PSU



**53506**  
0.6 inch x 0.2 inch  
(15 mm x 5 mm)  
PSU



**53510**  
0.6 inch x 0.32 inch  
(15 mm x 8 mm)  
PSU




**53518**  
0.6 inch x 0.47 inch  
(15 mm x 12 mm)  
PSU


Part#	Measurement	Material	Part#	Measurement	Material
53501	0.6 inch x 0.1 inch (15 mm x 2.5 mm)	PSU	53509	0.6 inch x 0.28 inch (15 mm x 7 mm)	PSU
53502	0.6 inch x 0.12 inch (15 mm x 3 mm)	PSU	53510	0.6 inch x 0.32 inch (15 mm x 8 mm)	PSU
53503	0.6 inch x 0.14 inch (15 mm x 3.5 mm)	PSU	53511	0.6 inch x 0.34 inch (15 mm x 8.5 mm)	PSU
53504	0.6 inch x 0.16 inch (15 mm x 4 mm)	PSU	53512	0.6 inch x 0.35 inch (15 mm x 9 mm)	PSU
53505	0.6 inch x 0.18 inch (15 mm x 4.5 mm)	PSU	53513	0.6 inch x 0.37 inch (15 mm x 9.5 mm)	PSU
53506	0.6 inch x 0.2 inch (15 mm x 5 mm)	PSU	53514	0.6 inch x 0.39 inch (15 mm x 10 mm)	PSU
53507	0.6 inch x 0.22 inch (15 mm x 5.5 mm)	PSU	53516	0.6 inch x 0.45 inch (15 mm x 11.5 mm)	PSU
53508	0.6 inch x 0.24 inch (15 mm x 6 mm)	PSU	53518	0.6 inch x 0.47 inch (15 mm x 12 mm)	PSU




**31300**  
0.6 inch x 0.1 inch  
(15 mm x 2.5 mm)  
PSU



**31306**  
0.6 inch x 0.22 inch  
(15 mm x 5.5 mm)  
PSU



**31310**  
0.6 inch x 0.3 inch  
(15 mm x 7.5 mm)  
PSU



**31316**  
0.6 inch x 0.45 inch  
(15 mm x 11.5 mm)  
PSU

Part#	Measurement	Material	Part#	Measurement	Material
31300	0.6 inch x 0.1 inch (15 mm x 2.5 mm)	PSU	31309	0.6 inch x 0.28 inch (15 mm x 7 mm)	PSU
31301	0.6 inch x 0.12 inch (15 mm x 3 mm)	PSU	31310	0.6 inch x 0.3 inch (15 mm x 7.5 mm)	PSU
31302	0.6 inch x 0.14 inch (15 mm x 3.5 mm)	PSU	31311	0.6 inch x 0.32 inch (15 mm x 8 mm)	PSU
31303	0.6 inch x 0.16 inch (15 mm x 4 mm)	PSU	31312	0.6 inch x 0.34 inch (15 mm x 8.5 mm)	PSU
31304	0.6 inch x 0.18 inch (15 mm x 4.5 mm)	PSU	31313	0.6 inch x 0.35 inch (15 mm x 9 mm)	PSU
31305	0.6 inch x 0.2 inch (15 mm x 5 mm)	PSU	31314	0.6 inch x 0.37 inch (15 mm x 9.5 mm)	PSU
31306	0.6 inch x 0.22 inch (15 mm x 5.5 mm)	PSU	31315	0.6 inch x 0.39 inch (15 mm x 10 mm)	PSU
31307	0.6 inch x 0.24 inch (15 mm x 6 mm)	PSU	31316	0.6 inch x 0.45 inch (15 mm x 11.5 mm)	PSU
31308	0.6 inch x 0.26 inch (15 mm x 6.5 mm)	PSU			



Luer Activated Inflation Valves



**80360**  
0.157 inch - 0.358 inch OD  
(4.0 mm - 9.1 mm)  
PP, PVC, Silicone



**80062**  
0.287 inch OD (7.3 mm)  
Nitrile, Polyester, PVC, SS



**51410**  
Assembly Cuff for 80063  
PP




**80063**  
0.181 inch OD (4.6 mm)  
Nitrile, Polyester, PP, SS



**80064**  
0.287 inch OD (7.3 mm)  
Nitrile, Polyester, PVC, SS



**12586**  
Bite Stick, White  
4.75 inch (121 mm)  
PE



**53523**  
Red Tracheal Plug with Tab  
Fits 6 FR OD Tubing (0.079 inch; 2 mm);  
0.073 inch ID (1.9 mm)  
Silicone



**56025**  
Double Swivel Elbow Connector  
without Suction Port  
15 mm ID x 24 mm OD  
15 mm ID  
K-Resin, PP



**51015**  
Double Swivel Elbow with Port/Plug  
15 mm ID x 15 mm OD  
PP, PU Blend and PVC



**56024**  
Double Swivel Elbow Connector with  
Suction Port  
15 mm ID x 23 mm OD  
15 mm ID  
K-Resin, PP, Thermoplastic Rubber



**13135**  
Straight Connector with Side Port  
15 mm OD x 15 mm ID/22 mm OD;  
Side Port Female Luer Lock  
PP



**53519**  
Side Port Autocontrol Airway Connector  
15 mm OD x 15 mm ID;  
Side Port 4.1 mm OD  
PP



**51306**  
Straight Adapter  
22 mm OD/15 mm ID x 22 mm ID;  
Side Port 4.1 mm ID  
PP



**53522**  
Disconnect Wedge  
15 mm  
PBT

Facilitates disconnection of tracheal tube  
from ventilator adapter

Minimizes tube torquing