



Spike Components and Caps

Spikes

Filter



**13125**  
0.256 inch ID (6.5 mm)  
3 Micron  
Microglass Fibers, PE



**X5008**  
0.291 inch Diameter (7.4 mm)  
3 Micron  
PP, PTFE



**X5009**  
0.217 inch Diameter (5.5 mm)  
3 Micron  
PP, PTFE



**23241**  
Air Vent  
3 Micron  
Fits 23226  
PP

Vented Spike Assemblies



**X1001**  
0.106 inch OD (2.7 mm)  
20 Drops/ml  
HDPE, PVC



Cap Included



**X1005**  
0.106 inch OD (2.7 mm)  
20 Drops/ml  
HDPE, PP, PTFE, PVC



Cap Included

Also Available:	
<b>X2002</b>	25 Micron

Also Available:	
<b>X2010</b>	25 Micron