

Female Luer Lock

23203

ABS, Acrylic Copolymer, PE, PP

25661 **NEW!**


ABS, Polyethylene, Versapor

11595

0.106 inch ID x 0.173 inch OD
(2.8 mm x 4.4 mm)
ABS, HDPE, PP

25660 **NEW!**

0.15 inch - 0.157 inch ID x
0.217 inch OD
(3.8 mm - 4 mm x 5.5 mm)
ABS, PE, Versapor



Cap Included

Cap Included

Cap Included

23205

0.161 inch ID x 0.189 inch OD
(4.1 mm x 4.8 mm)
ABS, Acrylic Copolymer, PE, PP

23202

0.15 inch - 0.169 inch ID
(3.8 mm - 4.3 mm)
ABS, Acrylic Copolymer, PE, PP

23209

0.209 inch OD (5.3 mm)
ABS, Acrylic Copolymer, HDPE, LDPE