



**NEWS RELEASE
FOR IMMEDIATE RELEASE**

Qosina's New Flexible Advanced Polymer Multifunctional Tie

Edgewood, N.Y., December 22, 2005 –Qosina is pleased to introduce [advanced polymer multifunctional ties](#). The tie design is based on modern cabling, a chain of cells terminated by a threading tongue. The tongue can be inserted through any cell, and the entire strip pulled through to form a closed loop. These releasable ties can be double wrapped or formed into multiple loops, adaptable to many different configurations. Additionally, the unused tail can be cut and reused as a new tie. These reusable, releasable ties are an economical, waste-free and versatile option for connecting tubing, cables or as a closure.

Qosina offers [Clear Flexible Ties \(#99912\)](#) and [UV Stable Black Flexible Ties \(#99913\)](#) in 12 inch lengths. Because of the unique structure of these ties, they can form closed loops of various sizes with multiple wraps in either direction. These very soft ties secure items with a light, even pressure. Thermoplastic polyurethane elastomer combined with the cell design helps cushion items from knocks and vibration. Unlike hook and loop strips, these ties are elastic enough to prevent over-tightening. Both Multifunctional Flexible Ties are available in packs of 100.

Visit qosina.com to see our broad selection of re-opening and non-reopening plastic ties, clamp, clips in many sizes, colors and materials. Call our Customer Service Specialists to obtain our 2005 catalog, request a sample or place an order at +1 (631) 242-3000, fax +1 (631) 242-3230 or e-mail info@qosina.com. Custom sourcing service is available through our knowledgeable specialists using Qosina's extensive resources.

Qosina is an ISO 9001:2000, ISO 14001 registered worldwide supplier of stock and custom-made OEM disposable components including: stock connectors, luers, caps, clamps, adapters, stopcocks, valves, spikes, syringes and tubing.

