



**NEWS RELEASE
FOR IMMEDIATE RELEASE**

Qosina Female to Female Luer Lock Connector

Edgewood, NY, May 12, 2006 –Qosina is pleased to add a new female to female luer lock (*Part # 13136*) to our inventory. The wide-winged connector securely fastens two male luer lock components together quickly and easily. This connector is composed of ABS which limits reagent reactions. Threads on the exterior collar provide torquing connections to the interior threads of male ended components. Assemble your kit or line without reconfiguring using our convenient double ended connector.

At Qosina, www.qosina.com, you can view our luers, luer locks, and connectors in addition to hose clamps, tubing and other components to fit all your design needs. Call our Customer Specialists to request our 2007 catalog, samples or to place an order at +1 (631) 242-3000, fax +1 (631) 242-3230 or e-mail info@qosina.com. Custom sourcing service is available through our knowledgeable Specialists using Qosina's extensive resources.

Qosina is an ISO 9001:2000, ISO 14001:2004 registered worldwide supplier of stock and custom-made OEM disposable components including: stock connectors, luers, caps, clamps, adapters, stopcocks, valves, spikes, syringes and tubing.

