



**NEWS RELEASE
FOR IMMEDIATE RELEASE**

Qosina Flow Regulator

Edgewood, NY, December 18, 2006 — Qosina announces a new economically priced flow regulator (**Part # 21292**). This white regulator has a distinguishable easy-to-turn dial with clearly marked blue flow adjustment indicators (0 to 250 ml/ hour). Our ABS flow regulator is 1-5/8 inches in length and is designed for use with 3mm ID x 4.1mm OD sized tubing. Qosina carries several different flow control components to meet your design needs.

Qosina's inventory contains tubing, luers, caps, spikes, bags, injection sites, spikes, clamps and other components to aid with your project requirements. Call our Customer Service Specialists to obtain our latest 2007 catalog, request a sample or place an order at +1 (631) 242-3000, fax +1 (631) 242-3230 or e-mail info@qosina.com. Visit us at www.qosina.com to see our latest products. Custom sourcing service is available through our knowledgeable specialists using Qosina's extensive resources.

Qosina is an ISO 9001:2000, ISO 14001:2004 registered worldwide supplier of stock and custom-made OEM disposable components including: stock connectors, luers, caps, clamps, adapters, stopcocks, valves, spikes, syringes and tubing.

